


























**Rhinoplasty  
Instruments  
Set**





Rhinoplasty Kit includes a number of extremely high quality instruments for rhinoplasty procedures so functional and cosmetics. The Golden Titanium Nitride coating indicates that these German Stainless Steel made precision instruments have been carefully designed to help improve the surgical result by reducing postoperative ecchymosis and edema.

Micro Osteotome with Rounded guard, <b>2.2mm 18cm</b>	
Micro Osteotome with Rounded guard, <b>3mm 18cm</b>	
Micro Osteotome with Rounded guard, <b>4mm 18cm</b>	
Micro Osteotome Double Edged Grinding <b>2mm 18cm</b>	
Micro Osteotome Double Edged Grinding <b>2.5mm 18cm</b>	
Micro Osteotome Double Edged Grinding <b>3mm 18cm</b>	
Micro Osteotome Double Edged Grinding <b>4mm 18cm</b>	
Osteotome Double Edged Grinding Rounded Corners, Recessed grip, Cut of Width <b>10mm 18cm</b>	
Osteotome Double Edged Grinding Rounded Corners, Recessed grip, Cut of Width <b>12mm 18cm</b>	
Osteotome Double Edged Grinding Rounded Corners, Recessed grip, Cut of Width <b>14mm 18cm</b>	
Osteotome Double Edged Grinding Rounded Corners, Recessed grip, Cut of Width <b>16mm 18cm</b>	
Osteotome Double Edged Grinding Rounded Corners, Recessed grip, Cut of Width <b>18mm 18cm</b>	

One prong Sharp Hook 15cm Length	
Two prong Sharp Hook 15cm Length	
Two Prong Hook, Sharp, 10mm Width 12cm	
Aufrecht Nasal Retractor, Width of blade 8.5mm Length of retractor blade 45mm, 15cm Length	
Glabella Rasp, Curved Push-Pull Cut Length 20cm	
Joseph Raspatory 20cm Length	
Bone Rasp, Tungsten Carbide Double Ended Rasp Blades <b>5 and 6 Medium</b> Length 21cm	
Bone Rasp, Tungsten Carbide Double Ended Rasp Blades <b>7 and 8 Medium</b> and <b>fine</b> Length 21cm	
Adson Brown Forceps 12cm Fine Side Grasping teeth	
Adson Dressing Forceps 12cm, Serrated	
Adson Tissue Forceps 1x2 Teeth 12cm	
Rhinoplastic Scissors, Tungsten Carbide Extra delicate straight Length 10cm	
Rhinoplastic Scissors, Tungsten Carbide Extra delicate Curved Length 10cm	

<p>Walter Angular Tungsten Carbide Scissors Blunt / Blunt 10cm Length ( Smooth End Scissors )</p>	
<p>Kilner Scissors Curved Flat End, Tungsten Carbide 14cm Length</p>	
<p><b>Tungsten Carbide</b> Needle Holder, 13cm Length</p>	
<p>Blakesly <b>Nasal Forceps</b>, Straight Size 1, Working Length = 11cm</p>	
<p>Fomon Dorsal <b>Nasal Scissors</b>, Supercut Angled, Serrated</p>	

<p>Giunta tungsten carbide +Supercut <b>Nasal scissors</b> curved 13.5cm</p>	
<p>Asch septum straightening forceps angled, 9" 23cm</p>	
<p>Sheen Cartilage Grid measuring, and photographing cartilage grafts. Measured in millimeters. Can be used as a cutting surface.</p>	
<p>Desmarres Lid Retractor Size 0, 5-1/2" (140mm) length</p>	
<p>Asch septum straightening forceps angled, 9" 23cm</p>	
<p>Cinelli guarded osteotome, 18.5CM 10MM, straight Stainless Steel. Double-guarded osteotome for hump removal.</p>	
<p>Silver Guarded Lateral Osteotome 18cm 1 Pcs Straight</p>	

<p>Silver Guarded Lateral Osteotome 18cm 1 Pcs Curved Right</p>	 <p>The image shows a silver guarded lateral osteotome, which is a surgical instrument with a long, curved blade and a textured handle. Below the photograph is a technical line drawing of the instrument, showing its profile and the guard mechanism.</p>
<p>Silver Guarded Lateral Osteotome 18cm 1 Pcs Curved Left</p>	 <p>The image shows a silver guarded lateral osteotome, which is a surgical instrument with a long, curved blade and a textured handle. Below the photograph is a technical line drawing of the instrument, showing its profile and the guard mechanism.</p>
<p>Ballenger Swivel Knife 7-1/2" 19 cm, 4 mm Blade</p>	 <p>The image shows a Ballenger Swivel Knife, which is a surgical instrument with a long, curved blade and a handle. The blade is shown in a swiveling position, and the handle has a textured grip.</p>
<p>Stainless Steel Sterilizing Box ! This set comes packed in this box</p>	 <p>The image shows a stainless steel sterilizing box, which is a rectangular container with a perforated lid and a handle. The box is shown open, revealing the interior.</p>